



# VIP<sup>®</sup>-300 Pupillometer

VIP



NEUR<sup>•••</sup>PTICS



## Very Accurately Measures Pupil Size under Various Light Conditions

Enables clinician to make best possible medical decision for LASIK and premium IOL patients based on pupil size, stated daily activities, and patient preferences.

- Portable, battery operated, hand-held device
- Very accurately measures pupil diameter ( $\pm 0.03$  mm)
- “**Dark**” setting, 2-second measurement
- “**Variable**” setting measures scotopic, low mesopic, and high mesopic pupil sizes in one 12-second sequence to simulate light conditions patient may experience in daily life
- Durable, ergonomic design
- Infrared camera, high precision optics, processor and LED light source



“The NeurOptics® Pupillometer...is very **compact, simple to use,** reasonably priced and **gives accurate automatic readings to 0.1 mm.** An accurate instrument for measuring the scotopic pupil is essential for refractive surgeons.”

—James Salz, M.D.

# VIP Technology Produces Reproducible Measurement Results

Testing has shown that the NeurOptics® Pupillometer is consistent from unit to unit and operator to operator. In fact, the NeurOptics Pupillometer has the highest accuracy and lowest error of all commercially available pupillometers<sup>1</sup>, while at the same time being the most economical hand-held infrared device.

The key to the reproducibility in measurement is the use of VIP® (Vertex Invariant Pupillometry®) technology developed by NeurOptics. This system compensates for the approximate 12 mm range in vertex distance<sup>2</sup> that may lead to inaccurate pupil measurement (error up to 20%), producing accurate measurements that are not affected by vertex distance<sup>3</sup>.

## Simple to Use and Operate

- One-button activation
- No calibration by the user
- Data is stored on the device and can be recalled and/or printed
- Video of the measurement can be played on the device's screen



### References:

<sup>1</sup>Schallenberg M, Bangre V, Steuhl KP, Kremmer S, Selbach JM., Comparison of the Colvard, Procyon, and Neuroptics pupillometers for measuring pupil diameter under low ambient illumination. J Refract Surg. 2010 Feb;26(2):134-43.

<sup>2</sup>Distance from the front surface of the eye to the back of a lens or optical device.

<sup>3</sup>The pupil is never at rest, and therefore a single, static, measurement is unlikely to be an accurate assessment of a pupil's true amplitude or range of diameter. To compensate for pupillary unrest, NeurOptics employs a dynamic measurement system which captures 30 pupil positions over an approximate 2-second scanning period, thus producing the weighted average pupil size and standard deviation.

## Ordering Information

| NeurOptics® VIP®-300  | Part Number     |
|---|-----------------|
| <b>System Includes:</b> VIP®-300 Pupillometer, Charging Station & Power Supply, Eye Cups (2), Carrying Case | VIP3-KIT-01     |
| Optional Accessories  | Part Number     |
| Printer Kit   | NEUR-PRTS445-BT |
| Barcode Scanner   | BCS-CC-01       |

**Caution:** Federal (USA) law restricts this device to sale by or on order of a physician. Refer to product package insert for instructions, warnings, precautions and complications.

## Bluetooth® Broadcast Range and Frequency

| Broadcast Function                                     | Range                                    | Frequency |
|--|--|-----------|
| Bluetooth Barcode Scanner to/from VIP-300 Pupillometer | Up to 100 yards depending on environment | 2.45 GHz  |

## Technical Specifications

| Parameter  | Description   |
|--|---|
| Measurement Characteristics  | Input = Human pupil sizing varying from 1 mm–9 mm                                   |
|  | Mean and standard deviation of pupil diameter at different background illuminations |
|  | Accuracy: +/- 0.03 mm   |
| Degree of protection against electric shock  | Pupillometer is double insulated (Class II protection)                              |
| Classification of the equipment against ingress of liquids   | Ordinary equipment  |
| Degree of safety of application in the presence of flammable anesthetic mixture with air or with oxygen or nitrous oxide | The equipment is not an AP or APG category equipment                                |
| Mode of Operation  | On Demand battery operation   |
| Power Supply   | Input: 100-240 VAC +/- 8%   |
|  | Output: 6V, 2.8 Amps  |
| Battery  | 3.7V 3350 mAmp/hour Li: Ion Cell  |
| Operating Environment  | Temperature Range: 18° C (65 F) to 30° C (86° F)                                    |
|  | Relative Humidity: 20% to 70% RH. Non condensing at all times                       |
| Transportation and storage environment   | Temperature Range: 0° C (32° F) to 75° C (167° F)                                   |
|  | Relative Humidity: 10% to 95% RH. Non-condensing at all times                       |
| Dimensions   | With eye cup = 7.5" H, 3.5" W, 4.5" D   |
|  | Without eye cup = 7.5" H, 3.5" W, 3.5" D  |
| Weight   | 320 grams +/- 10 grams  |
| Classification   | Class 1 LED product per IEC 60825   |

© 2016 NeurOptics, Inc. NeurOptics and VIP are all trademarks of NeurOptics, Inc. Bluetooth is a registered trademark of Bluetooth SIG, Inc.



18101 Von Karman Avenue, Suite 1940  
Irvine, CA 92612 | USA  
p: 949.250.9792  
Toll Free North America: 866.99.PUPIL  
info@NeurOptics.com  
NeurOptics.com

MKG-VIP3-01 Rev A